



High Resolution Epitope Mapping and Antibody Cross-Reactivity Analysis

Carsten Haber¹, Volker Stadler¹


¹ PEPperPRINT GmbH, Rischerstraße 12, 69123 Heidelberg, Germany

PEPPERPRINT
A NEW DIVERSITY

Abstract

Antibodies are among the most important life science tools for therapeutics, basic research and diagnostic tests. However, mono- and polyclonal antibodies are often poorly characterized in terms of specificity and cross-reactivity, needing further validation and cross-reactivity analysis. High density peptide microarrays are ideally suited for antibody characterization, enabling the analysis of epitope-antibody interactions with unmet speed and precision. Based on approaches for linear and conformational epitope mappings as well as high resolution epitope substitution scans, we developed a comprehensive toolbox for the in-depth analysis of epitopes and antibody cross-reactions for all kinds of antibodies and isotypes in the most flexible, comprehensive and economic manner. Moreover, we designed a Human Epitome Microarray comprising all linear human B-cell epitopes of the Immune Epitope Database (23,163 linear peptides), complemented by 4,661 epitopes of the most common vaccines. In a three-step approach based on a single assay, the Human Epitome Microarray enables a very detailed cross-reactivity analysis of antibodies including the identification of antibody-specific consensus motifs and database blasting to identify possibly cross-reacting human antigens.

Characterization of Rituximab by Conformational Epitope Mapping and Epitope Substitution Scan



- Rituximab (Rituxan®), a chimeric monoclonal anti-CD20 antibody, is known to interact with the 15 amino acid loop of the CD20 antigen (Du et al., 2007): NIYNCEPANPSEKNSPSTQYCYSIQ (Figure 1).
- High resolution epitope mapping with peptide microarrays carrying linear and cyclic constrained overlapping peptides (7 aa, 10 aa and 13aa) of the CD20 antigen.

Figure 1: Structure of the Rituximab Fab-CD20 epitope-peptide complex (Du et al., 2007).

Results Epitope Mapping

- Linear CD20 peptides did not show any response with Rituximab (data not shown), however constrained cyclic CD20 peptides show a very clear and strong response.
- High resolution data identifying **EPANPSEK as Rituximab epitope**. Clear epitope peaks for constrained cyclic peptides with excellent signal to noise ratio (Figure 2).




Figure 2: Peptide Microarray of cyclic CD20 peptides. Array was incubated with Rituximab followed by staining with anti-human IgG (red) and anti-Flag (green) antibodies. First double-row in each half corresponds to the 7 amino acid peptides, the second and the third double-row to the 10 and 13 amino acid peptides. Intensity plots of the Rituximab assays with linear and constrained cyclic CD20 peptides (right).

Results Epitope Substitution Scan

- In-depth epitope analysis of the constrained cyclic CD20 peptide NIYNCEPANPSEK led to a clear and well-defined substitution pattern.
- Array and amino acid plot reflect the conserved core motif EPANPSEK as well as a variable amino acid stretch NIYNC at the N-terminus (Figure 3).




Figure 3: Epitope Substitution Scan of wild type Rituximab peptide NIYNCEPANPSEK (left). Each amino acid position was gradually exchanged by the 20 L-amino acids. The same substitution pattern is visualized by an amino acid plot (right).

Antibody Cross-Reactivity Profiling of a Human Monoclonal Antibody with the PEPperCHIP® Human Epitome Microarray

- Cross-reactivity analysis of a human monoclonal anti-c-Myc antibody with the PEPperCHIP® Human Epitome Microarray in a three-step approach (Figure 4):
 - The Human Epitome Microarray covers 29,127 linear peptides printed in duplicate (58,254 peptide spots) as well as additional polio and HA control peptides.
 - Using the MEME tool (<http://meme-suite.org/tools/meme>), common motifs in the peptide top hits can be discovered.
 - Subsequent FIMO analysis (<http://meme-suite.org/tools/fimo>), translates hit motifs into possible cross-reactive antigens.





Figure 4: Workflow of PEPperPRINT's three-step approach for antibody cross-reactivity analysis.


Results Antibody Cross-Reactivity Profiling

- PEPperCHIP® Human Epitome Microarray incubated with human mAb followed by staining with secondary and control antibodies resulted in few but clear peptide hits (red).
- Cross-reactions against epitopes of nucleoprotein, erythrocyte membrane-associated giant protein antigen 332 and dihydrolipoamide S-acetyltransferase.
- Cross-reactivity profile over 28,895 known human epitopes with database annotations.
- Identification of conserved and variable amino acid positions by MEME motif discovery.

Figure 5: PEPperCHIP® Human Epitome Microarray of human mAb anti-c-Myc. Interactions of anti-c-Myc with the database epitopes in red and well-defined frame of HA control peptides in green (left). Top 40 interactions of human mAb anti-c-Myc sorted by decreasing spot intensities. MEME analysis of the top 10 peptide hits of the PEPperCHIP® Human Epitome Microarray (right).




Technology



- High-density PEPperCHIP® peptide microarrays are generated by digital laser printing on standard glass slides using a proprietary laser printer comprising 24 cartridges filled with individual amino acid toners.
- For array production, the amino acid toner are simultaneously printed with high precision on their respective positions on the glass slides. Peptide synthesis is initiated by melting the toner particles to 90 °C. Under these conditions, the amino acids are released and are available for coupling to the previous amino acid. The coupling cycle is completed by washing steps to remove excess building blocks and protecting groups. Finally, the array is ready for the next synthesis cycle with laser printing and coupling.
- The benefits of this technique are an unique flexibility in terms of peptide content, a high spot density with up to 11,000 features per chip and low material consumption enabling the generation of customized peptide array at reasonable costs.


PEPperCHIP® Peptide Microarray Platform



Antibody-peptide interactions are analyzed by immuno-type assays in a high-throughput fashion on peptide microarrays. Depending on the application, various microarray formats are available:

- Discovery format:** approximately 70,000 individual peptides; suitable for screening of large, diverse epitope libraries or several full-length viral proteomes on a single chip; applied for epitope, biomarker and target binder discovery as well as antibody cross-reactivity analysis
- Standard format:** covers approximately 11,000 individual peptides; routinely used for epitope mapping. Custom peptide microarrays or standard chips such as PEPperCHIP® Infectious Disease Epitope Microarray or the PEPperCHIP® Dengue Virus Type 1 Proteome Microarray among others
- Mapping format:** several identical array copies on a single chip; ideal for parallel screening of multiple samples; used for epitope mapping of single protein antigens, detailed epitope characterization or biomarker validation
- Multiplexed format:** up to 16 array copies on a single chip; ideal for assay development or hit validation studies with sample cohorts

Summary and Outlook



- High Resolution Epitope Mapping of Rituximab** with cyclic constrained peptides reveals a conserved core motif EPANPSEK
- Cyclic constrained peptides on the peptide microarray enable the specific detection of the Rituximab binding signal.
- Loop sizes of cyclic constrained peptides influence binding signal intensity.
- Characterization of other, unknown epitopes of therapeutic and diagnostic antibodies can be conducted with high density peptide microarrays.
- Antibody characterization and cross-reactivity analysis** of a human monoclonal antibody (anti-c-Myc (ch19E10)) is possible with the PEPperCHIP® Human Epitome Microarray - based three step approach.
- Antibody Cross-Reactivity Profiling** reveals possible cross-reactivity with antigens containing the hit motif 'XLV(S/A/P)E'. Top candidates were based on the consensus motif HLVS(A/P)E and assigned to DNA repair protein XRCC4, putative uncharacterized protein BVES-AS1, transmembrane protein 109 or choline O-acetyltransferase.
- Identification of other therapeutic and diagnostic antibodies as well as cross-reactivity analysis can be performed with PEPperCHIP® Human Epitome Microarray.