

Immunological Services | Regulatory T-Cells

Analyze Treg-mediated immune responses to validate your therapeutic candidates

Regulatory T-cells (Tregs) are essential for controlling the immune response and maintaining self-tolerance. Hence, **Tregs are promising cellular targets for cancer and autoimmune diseases treatment**. While preventing Treg-mediated immunosuppression can help to reinforce antitumor immune responses, forcing an induced Treg (iTregs) phenotype can temper pathologies in autoimmune diseases or transplant rejection.

With our **iTreg phenotyping** service, we can monitor the success of an ongoing clinical trial. Furthermore, via our **iTreg polarization** assay and **iTreg suppression** assay, we can support the development of Treg-targeting drug candidates by validating the modulating effect of your therapeutic candidates on T-cell polarization and function. The iTreg polarization and iTreg suppression assays are established for drug screenings in immune-oncology and autoimmune diseases.

Applications


Monitor *in vitro* the amount and character of Tregs during **patients' treatment in clinical trials**.

Test the effectiveness of a drug in enhancing the regulatory phenotype and functional character of iTregs, which could be beneficial for **treating autoimmune diseases or transplant rejection**.

Revise the efficiency of a drug candidate by preventing Treg differentiation and inhibiting iTreg function. This can support cancer treatment by inhibiting Tregs to boost antitumor immune responses.


Treg Phenotyping Assay Service

The Treg phenotyping assay is based on antibody-staining of Treg-specific transcription factors, surface markers, and quantitative analyses by flow cytometry. It allows the **characterization of Tregs** among the T-lymphocyte population or peripheral blood mononuclear cell (PBMC) samples. Subsequently, the identified Treg population can be sorted for further functional assays.




iTreg Polarization Assay Service

The iTreg polarization assay is a rapid and reliable *in vitro* system to quantitatively validate the capability of therapeutic candidates to modulate differentiation of naive T lymphocytes towards an iTreg phenotype, or to verify the candidate compound-mediated prevention of iTreg polarization.




iTreg Polarization Assay. Naive CD4⁺ T cells isolated from peripheral blood mononuclear cells (PBMC) are activated with a Treg-polarization cytokine cocktail. The drug candidates can be screened for their ability to boost or prevent the differentiation of Tregs via analyzing the expression of Treg-specific markers such as FOXP3, CD127, and CD25 via flow cytometry.



Representative FACS plots of a differentiation assay showing successful generation of iTregs from naive CD4⁺ cells. FOXP3 and CD25 expression in control T-cells and iTregs induced by the differentiation protocol above is shown. The pseudocolor plots show representative FOXP3 and CD25 stainings for one donor, pre-gated on singlet, live CD4⁺ cells.

iTreg Suppression Assay Service

The iTreg suppression assay constitutes a rapid and reliable *in vitro* functional test that validates the iTreg-induced suppression of functional T-cells by analyzing the proliferation of cocultured T lymphocytes. The assay can **monitor the modulating effect of drug candidates on Treg induction** and reveal booster or suppressor functions. This assay complements the iTreg polarization assay to validate the T-cell modulating effect of the drug candidate.



iTreg Suppression Assay. Fluorescently labeled T effector cells are activated and cultured alone or with polarized autologous iTreg cells and treated with the drug candidates. The proliferation of T effector cells and the suppressive effect of iTregs on the proliferation is analyzed by measurement of the decreasing fluorescence signal via flow cytometry. Finally, the modulation of iTreg-mediated suppression through the drug candidate is assessed.