

## Application Note

# Discovery of Overlapping Immunogenic B-cell and T-cell Epitopes

### ABSTRACT

This application note describes the use of PEPperCHIP® Peptide Microarrays combined with functional T-cell activation assays to discover coincident B-cell receptor and T-cell receptor epitopes (Figure 1). We applied this approach to unravel immunogenic epitopes in Epstein-Barr virus nuclear antigen-1 (EBNA-1) that activate B-cells to produce antibodies as well as T-cells to produce interferon gamma (IFN- $\gamma$ ). For this purpose, we analyzed sera and peripheral blood mononuclear cells (PBMCs) from two different donors. High resolution epitope mapping with the PEPperCHIP® Epstein-Barr Virus Peptide Microarray enabled the identification of promising B-cell receptor epitopes in EBNA-1 that were further evaluated in terms of T-cell activation in an IFN- $\gamma$  ELISpot assay. This approach allowed the discovery of a peptide sequence that functions as B-cell as well as T-cell immunogenic epitope.

### INTRODUCTION

Due to the huge B-cell and T-cell receptor repertoire, the lymphocyte-driven humoral and cellular immune response is highly specific (1). Although a foreign protein contains a large variety of different peptides that could evoke an immune response, only a few peptides are actually immunogenic and even less are immunodominant (2). Several peptide characteristics such as peptide length, amino acid sequence and human leukocyte antigen (HLA) affinity constitute a peptide as lymphocyte receptor epitope. Importantly, a B-cell epitope might also function as T-cell epitope. However, for most peptide sequences this is not the case (2).


Depending on the research topic, the simultaneous activation of B-cells and T-cells can be desired. In the development of a protective vaccines that induce contemporary as well as persisting immunity through humoral and cellular responses, overlapping epitopes for B-cell and T-cell receptor would be highly favorable.

For the selection of immunogenic peptide candidates, a variety of epitope prediction tools (such as the prediction tools from IEDB Analysis Resource) can help to forecast immunogenicity by using published data and/or applying

knowledge on preferred amino acids lengths and sequences etc. (3). However, the veracity of predictions is unambiguously limited and cannot replace laboratory screenings of B-cell and T-cell receptor epitopes.

High-density peptide microarrays are a powerful tool to simultaneously screen tens of thousands of peptides against serum antibodies in a high-throughput context (4). Subsequently, identified B-cell receptor epitopes can serve as basis for the generation of synthetic overlapping peptide libraries for T-cell immunogenicity testing by ELISpot assays. ELISpots directly determine the frequency of antigen-specific T-cells by detecting released activation markers such as IFN- $\gamma$ . The method is extremely sensitive and enables a quantitative single cell analysis that is impossible with supernatant-based assays such as ELISAs. Moreover, multiparametric ELISpot (or FluoroSpot) assays allow the detection of different markers in parallel and thereby the characterization of T-cell subtypes.

This application note describes an experimental workflow for exploring coincident B-cell receptor and T-cell receptor epitopes in EBNA-1 protein of EBV. B-cell receptor epitopes were identified with PEPperCHIP® Epstein-Barr Virus Peptide Microarrays. The corresponding peptides were incorporated in an overlapping mini peptide library and



**Figure 1. Workflow for the discovery of coincident B-cell and T-cell epitopes.** Serum samples are screened with PEPperCHIP® Peptide Microarrays displaying peptide libraries of up to 35,000 different peptides printed in duplicate. Identified hit peptides are validated via literature review and T-cell activation prediction tools. Subsequently, promising peptide candidates are synthesized applied in ELISpot/FluoroSpot assays with PBMCs of the same donor.

tested for T-cell activation via IFN- $\gamma$  ELISpot. Via this approach we were able to identify FHPVGEADY in EBNA-1 as a core sequence of B-cell as well as T-cell receptor epitope.

## RESULTS AND DISCUSSION

In the study, we initially mapped the humoral immune responses of EBV-infected individuals. To determine the infection-elicited antibodies against EBNA-1, we screened sera of two donors with PEPperCHIP® Epstein-Barr Virus Peptide Microarrays (Figure 2). The PEPperCHIP® Epstein-Barr Virus Peptide Microarrays contains 5,549 linear peptides of the most immunogenic EBV antigens printed in duplicates, and were framed by poliovirus VP1 protein 20-22 and influenza hemagglutinin (HA) derived control peptides. The microarray also includes 257 overlapping EBNA-1 peptides for high-resolution epitope data.


PEPperCHIP® Epstein-Barr Virus Peptide Microarrays were initially incubated with the secondary goat anti-human IgG (Fc) DyLight680 antibody to identify possible background interactions of the detection antibodies. Here, no background interactions were observed. Subsequently, the sera of donor 1 and donor 2 were diluted to 1:150 and incubated overnight at 4°C. After washing, detection was done using the secondary goat anti-human IgG (Fc) DyLight680

antibody. Fluorescence readout was performed using an INNOPSYS Imaging System, and 16-bit gray-scale images were analyzed with the PepSlide® Analyzer software.

The in-depth analysis of the IgG responses highlighted clearly differential antibody profiles of both donors against EBNA-1. One of the main epitopes of donor 1 was based on peptides with the consensus motif FHPVGEADY (individual peptides RRPFFHPVGEADY, PFFHPVGEADYFE, FHPVGEADYFEYH). In contrast, the serum of donor 2 did not show remarkable antibody responses against the EBNA-1 peptides and particularly with epitope FHPVGEADY (Figure 2).

To discover potential coincident T-cell receptor epitopes within the identified B-cell receptor epitopes, peptides RRPFFHPVGEADY, PFFHPVGEADYFE and FHPVGEADYFEYH were synthesized. In addition, peptides HPVGEADYFEY of EBNA-1 protein and DYC NVLNKEF from replication and transcription activator (BRLF1) were included in the study. Peptide HPVGEADYFEY is described to be highly immunogenic in EBV-exposed individuals expressing HLA-B\*3501 (5-7), while peptide DYC NVLNKEF is known to activate T-cells in HLA-A\*2402-expressing individuals (8).

Peptide-mediated T-cell activation was measured via interferon gamma (IFN- $\gamma$ ) ELISpot. IFN- $\gamma$  is commonly used as T-cell activation marker in ELISpot and has a wide range of functions in the immune response, such as



**Figure 2. Discovery of B-cell receptor epitopes via peptide microarray. (A) Scans of the PEPperCHIP® Epstein-Barr Virus Peptide Microarrays and (B) fluorescence intensity plots.** The microarrays were incubated with the sera of two donors at a dilution of 1:150 overnight at 4°C. Detection was done using the secondary goat anti-human IgG (Fc) DyLight680 antibody. Fluorescence readout was performed using an INNOSYS Imaging System. Red spots=IgG responses. The frame corresponds highlights the position of the EBNA-1 peptides. **(C) Antibody response profile against EBNA-1.** Scanned images were analyzed with the PepSlide® Analyzer software. The intensity plots show the fluorescence intensities of 257 overlapping EBNA-1 peptides sorted from the N- to the C-terminus of the protein. While we observed moderate IgG responses of donor 1 with two clear epitopes, the serum of donor 2 exhibited only weak signals and no clear epitope.

activation of macrophages and induction of IgG production from activated B cells. IFN- $\gamma$  is mainly produced and secreted by activated CD4 and CD8 T-cells which constitute about 45–70% of PBMCs in human peripheral blood (9,10).

Briefly, PBMCs were tested directly in a 96-well IFN- $\gamma$  ELISpot plate for 24 hours. 500,000 PBMCs per well were stimulated in triplicates either with the mini peptide library made of the 5 different synthetic (Table 1) or with the single peptides in a final concentration of 10  $\mu$ g/ml. As positive control, PBMCs were activated with a commercial EBNA-1 peptide mix containing 158 15 amino acid peptides with 11 amino acids peptide-peptide overlap (JPT, Berlin, Germany).


The assay was performed according to the manufacturer's protocol and analyzed on an ImmunoSpot® Analyzer (CTL ImmunoSpots Europe, Bonn, Germany). As expected, the positive control EBNA-1 peptide mix induced a pronounced IFN- $\gamma$  response in both donors. Donor 1 could also be activated with the mini peptide library. Based on the tests with the single peptides, this activation was due to immunogenicity of peptides FHPVGEADYFEYH, HPVGEADYFEY and DYC NVLNKEF. To a minor degree, donor 2 could also be activated with the mini peptide library. Interestingly, this

activation was exclusively due to the immunogenicity of peptide DYC NVLNKEF, as none of the EBNA-1 derived peptides induced an IFN- $\gamma$  response (Figure 3). This result was in line with the peptide microarray data that highlighted a humoral response against FHPVGEADY-containing peptides in donor 1, but no corresponding response in donor 2. EBNA-1 peptides FHPVGEADYFEYH and HPVGEADYFEY induced T-cell activation in donor 1, whereas EBNA-1 peptides RRPFFHPVGEADY and PFFHPVGEADYFE did not. Therefore, the core T-cell epitope seemed to depend on appropriate C-terminal sequences.

**Table 1. Peptide sequences for validation of T-cell immunogenicity.**

Peptide sequence	Protein
RRPFFHPVGEADY PFFHPVGEADYFE FHPVGEADYFEYH HPVGEADYFEY	EBNA-1
DYC NVLNKEF	BRLF1

Peptides RRPFFHPVGEADY, PFFHPVGEADYFE and FHPVGEADYFEYH were identified as B-cell epitopes based on the previous epitope mapping with the serum samples of donor 1 and 2. Peptides HPVGEADYFEY of EBNA-1 protein and DYC NVLNKEF of replication and transcription activator (BRLF1) were chosen from literature. For T-cell stimulation, all five peptides were either applied simultaneously as mini peptide library or as individual peptides.



**Figure 3. Discovery of T-cell receptor epitopes via functional assay.** 500,000 PBMCs per well were stimulated in a 96-well IFN- $\gamma$  ELISpot plate for 24 hours with a commercial EBNA-1 peptide mix (EBNA-1), the mini peptide library (Table 1) or the underlying single peptides in a final concentration of 10  $\mu$ g/ml. ELISpot was analyzed on an ImmunoSpot<sup>®</sup> Analyzers. Top half: Quantification of ELISpots with the mean values and the standard deviation of triplicates of the counted spots per 1 million PBMCs. Bottom half: IFN- $\gamma$ -dependent blue-colored immune complexes.

## CONCLUSION

In conclusion, we were able to demonstrate that the combination of PEPperCHIP<sup>®</sup> Peptide Microarrays for antibody analysis and functional T-cell activation assays is ideally suited to unravel coincident B-cell receptor and T-cell receptor epitopes. The approach can be easily transferred to any other immunogenicity testing with custom antigens e.g. for vaccine development, cancer research or the analysis of anti-drug responses.

## BIBLIOGRAPHY

- Jackson KJ, Kidd MJ, Wang Y, Collins AM. The shape of the lymphocyte receptor repertoire: lessons from the B cell receptor. *Front Immunol.* 2013 Sep 2;4:263. doi: 10.3389/fimmu.2013.00263. PMID: 24032032; PMCID: PMC3759170.
- Malonis RJ, Lai JR, Vergnolle O. Peptide-Based Vaccines: Current Progress and Future Challenges. *Chem Rev.* 2020 Mar 25;120(6):3210-3229. doi: 10.1021/acs.chemrev.9b00472. Epub 2019 Dec 5. PMID: 31804810; PMCID: PMC7094793.
- Banerjee S, Majumder K, Gutierrez GJ, Gupta D, Mittal B. Immuno-informatics approach for multi-epitope vaccine designing against SARS-CoV-2. *bioRxiv [Preprint].* 2020 Aug 17:2020.07.23.218529. doi: 10.1101/2020.07.23.218529.
- Heiss K, Heidepriem J, Fischer N, Weber LK, Dahlke C, Jaenisch T, Loeffler FF. Rapid Response to Pandemic Threats: Immunogenic Epitope Detection of Pandemic Pathogens for Diagnostics and Vaccine Development Using Peptide Microarrays. *J Proteome Res.* 2020 Nov 6;19(11):4339-4354. doi: 10.1021/acs.jproteome.0c00484.
- Miles, J.J., N.A. Borg, R.M. Brennan, F.E. Tynan, L. Kjer-Nielsen, S.L. Silins, M.J. Bell, J.M. Burrows, J. McCluskey, J. Rossjohn, and S.R. Burrows. 2006. TCR alpha genes direct MHC restriction in the potent human T cell response to a class I-bound viral epitope. *J. Immunol.* 177:6804–6814.
- Tellam, J., G. Connolly, K.J. Green, J.J. Miles, D.J. Moss, S.R. Burrows, and R. Khanna. 2004. Endogenous presentation of CD8+ T cell epitopes from Epstein-Barr virus–encoded nuclear antigen 1. *J. Exp. Med.* 199:1421–1431. doi: 10.1084/jem.20040191
- Lee, S.P., J.M. Brooks, H. Al-Jarrah, W.A. Thomas, T.A. Haigh, G.S. Taylor, S. Humme, A. Schepers, W. Hamerschmidt, J.L. Yates, et al. 2004. CD8 T cell recognition of endogenously expressed Epstein-Barr virus nuclear antigen 1. *J. Exp. Med.* 199:1409–1420. doi: 10.1084/jem.20040121
- Rafiee N, Ravanshad M, Asadi B, Kianfar R, Maleki A: Investigation of IL-2 and IFN- $\gamma$  to EBV Peptides in Stimulated Whole Blood among Multiple Sclerosis Patients and Healthy Individuals. *Intervirolgy* 2021;64: 203-208. doi: 10.1159/000517002

9. Slota M, Lim JB, Dang Y, Disis ML. ELISpot for measuring human immune responses to vaccines. *Expert Rev Vaccines*. 2011 Mar;10(3):299-306. doi: 10.1586/erv.10.169.
10. Yang F, Patton K, Kasprzyk T, Long B, Gupta S, Zoog SJ, Tracy K, Vettermann C. Validation of an IFN-gamma ELISpot assay to measure cellular immune responses against viral antigens in non-human primates. *Gene Ther*. 2022 Feb;29(1-2):41-54. doi: 10.1038/s41434-020-00214-w.

 PEPperPRINT GmbH  
Tullastr. 2  
69126 Heidelberg  
Germany

 [www.pepperprint.com](http://www.pepperprint.com)  
 [info@pepperprint.com](mailto:info@pepperprint.com)  
 +49-6221-7264489

 [pepperprint-gmbh](https://www.linkedin.com/company/pepperprint-gmbh)  
 [@pepperchip](https://twitter.com/pepperchip)  
 [PEPperPRINTGmbH](https://www.facebook.com/PEPperPRINTGmbH)  
 [PEPperPRINT](https://www.youtube.com/PEPperPRINT)

