



### Day 1



### Day 2



## Protocols

### Download and read

[pepperprint.com/technology/downloads-and-protocols](http://pepperprint.com/technology/downloads-and-protocols)

### Watch

- Animated tutorial: [youtu.be/-kc5MV8z8GY](https://youtu.be/-kc5MV8z8GY)
- JoVE video: [youtu.be/sY6vFbTBWoU](https://youtu.be/sY6vFbTBWoU)

## Buffers

### Wash buffer

PBS, pH 7.4  
+ 0.05% Tween 20 for linear peptides  
+ 0.005% Tween 20 for cyclic constrained peptides

### Blocking buffer

Rockland Blocking buffer MB-070  
or  
Wash buffer with 1% BSA

### Staining buffer

Wash buffer with 10% blocking buffer

### Dipping buffer

1 mM Tris, pH 7.4

## General recommendations

- Only touch the microarray slide at the rims with your fingers, or use plastic forceps (in the area of microarray number).
- Place the microarray number appearing mirror-inverted in the upper right corner.
- Use the dark incubation tray lid to prevent dye bleaching.
- Use an orbital shaker for optimal surface wetting (avoid rocking incubation, which might increase background and artefacts).
- Pipette with the tip being placed close to the incubation tray wall; avoid touching the microarray content.



Approximate experimental time for each day: 2.5-3 hours

Watch our video tutorial:



[youtu.be/-kc5MV8z8GY](https://youtu.be/-kc5MV8z8GY)

PEPperPRINT GmbH  
Tullastr. 2  
69126 Heidelberg  
Germany

[www.pepperprint.com](http://www.pepperprint.com)  
[info@pepperprint.com](mailto:info@pepperprint.com)  
+49-6221-7264489

### Sample volumes and dilutions

Sample type or secondary antibody	Range	Recommended starting dilution or concentration	PEPperCHIP® array format (no. of array copies/chip)	Minimum assay volume
Purified antibody	1 - 500 µg/mL	1 µg/mL	1	1,500 µL
Serum (IgG, IgM, IgA)	1:10 - 1:10,000	1:1,000	2	700 µL
Serum (IgE)	1:2 - 1:10,000	1:5	3	400 µL
Cerebrospinal fluid	1:10 - 1:10,000	1:100	4	350 µL
Secondary antibody	1:500 - 1:10,000	1:100	5	200 µL
Anti-HA control antibody	1:2,000	1:2,000	16	100 µL

### Troubleshooting

Issue	Possible reason	Solution
<b>High background signal</b>	Insufficient washing	Extend washing steps.
	Non-homogenous sample	Centrifuge sample prior to dilution.
	Use of strongly charged fluorescent dye	Switch to a combination of either Cy3/Cy5 or DL680/DL800.
<b>Weak signal intensities</b>	Low-affinity interactions or high off-target binding rates	Reduce washing steps or re-incubate the same array copy with a lower dilution (for serum) or a higher concentration (for antibody) of sample.
	Insufficient blocking	Use the recommended Rockland Blocking buffer for optimal signal-to-noise ratio.
<b>Artefacts on the microarray scan</b>	Portions of the microarray were left uncoated with sample or buffer	Ensure that the array is always fully coated with sample or buffer during the immunoassay.
	Array was poorly rinsed and dried	Dip the microarray chip according to the protocol and dry it immediately with a pressurized air stream.
	Incorrect placement of the silicone seal	Place and align the silicone seal precisely onto the microarray chip.
<b>Scratches or blurry areas on the microarray scan</b>	Accidental contact with the printed surface of the microarray chip	Avoid any contact to the microarray surface; only touch the microarray chip at the rim or use plastic forceps for handling.